

Skin Lesion Analysis

By Umut Celik Email: umut.celik@cix.csi.cuny.edu

CSC 741 - Digital Image Processing Midterm Project

1. Introduction

The Challenge

- Automated analysis of dermoscopic images is crucial for early skin cancer detection.
- These images present challenges: varying illumination, hair, diverse skin tones, and subtle lesion characteristics.

Project Goal

- To develop an image processing pipeline to:
 1. Segment (isolate) skin lesions from dermoscopic images.
 2. Extract a comprehensive set of visual features from the segmented lesions.
- Apply fundamental Digital Image Processing techniques to a real-world medical imaging problem.

Dataset

- **Source:** SIIM-ISIC Melanoma Classification (Kaggle).
 - **Content:** Dermoscopic images (JPEG format) and associated metadata (e.g., labels).
-

2. Image Processing Pipeline

A multi-step approach to process images and extract meaningful data:

Image Loading & Conversion -> Preprocessing -> Color Space Transformation -> Segmentation -> Feature Extraction

Step 1: Image Loading & Initial Conversion

- **Action:** Load JPEG images.
 - **Conversion:** Convert from BGR (OpenCV default) to RGB for consistency in processing and display.
 - Files involved: `src/data_loader.py`
-

Step 2: Preprocessing

- **Grayscale Conversion:**
 - Convert RGB images to grayscale for intensity-based operations.
 - Function: `rgb_to_grayscale` in `src/color_utils.py`
 - **Hair Removal:**
 - Technique: Morphological Black-Hat filtering to identify hair structures, followed by inpainting to remove them.
 - Reduces noise and artifacts caused by hair.
 - Function: `remove_hair` in `src/preprocessing.py`
 - **(Illumination Correction):**
 - Concept: Use morphological opening to estimate and correct non-uniform background illumination.
 - Status: Implemented (`correct_illumination` in `src/preprocessing.py`), but noted as potentially skipped in some display pipelines for directness to segmentation.
-

Step 3: Color Space Transformation

- **RGB to HSV:**
 - Convert RGB images to Hue, Saturation, Value (HSV) color space.
 - HSV is often more intuitive for color-based feature extraction:
 - **H (Hue):** Dominant color.
 - **S (Saturation):** Purity/intensity of color.
 - **V (Value):** Brightness.
 - Function: `rgb_to_hsv` in `src/color_utils.py`
-

Step 4: Lesion Segmentation

- **Goal:** Create a binary mask that accurately isolates the lesion area.
 - **Method: Otsu's Thresholding**
 - Applied to the preprocessed (hair-removed) grayscale image.
 - Automatically determines an optimal threshold value to separate lesion pixels from background.
 - *Function: `otsu_threshold` in `src/segmentation.py`*
 - **Mask Refinement:**
 - Technique: Morphological operations.
 - **Opening:** Erosion followed by dilation (removes small noise/objects).
 - **Closing:** Dilation followed by erosion (fills small holes within the lesion).
 - Improves the quality and contiguity of the mask.
 - *Functions: `opening`, `closing` in `src/morphology.py`, utilized by `apply_threshold` in `src/segmentation.py`*
-

Step 5: Feature Extraction

Extracted from the segmented lesion area (defined by the refined mask):

Feature Extraction: Intensity Features

- **From Grayscale image:**
 - Mean pixel intensity
 - Standard deviation of pixel intensities
 - Minimum pixel intensity
 - Maximum pixel intensity
 - *Module: `calculate_intensity_stats` in `src/feature_extraction.py`*
-

Feature Extraction: Color Features (HSV)

- **From HSV channels:**
 - **Hue (H):**
 - Mean Hue (using circular mean calculation for angular data)
 - Standard deviation of Hue
 - **Saturation (S):**
 - Mean Saturation
 - Standard deviation of Saturation
 - **Value (V):**
 - Mean Value (Brightness)
 - Standard deviation of Value
 - *Module: `calculate_hsv_stats` in `src/feature_extraction.py`*
-

Feature Extraction: Border Features

- **From Grayscale image & mask:**
 - Calculated using morphological gradient on the lesion border.
 - Features include:
 - Mean, Standard Deviation, and Max of border gradient values.
 - Border Irregularity (ratio of border pixels to mask area).
 - *Module: `calculate_gradient_features` in `src/border_texture_utils.py`*
-

Feature Extraction: Texture Features

- **From Grayscale image & mask:**
 - Calculated using Morphological Top-Hat (highlights bright details) and Bottom-Hat (highlights dark details) transforms.
 - Features include:
 - Mean and Standard Deviation of Top-Hat transform values.
 - Mean and Standard Deviation of Bottom-Hat transform values.
 - Texture Contrast Index (sum of Top-Hat and Bottom-Hat means).
 - *Module: `calculate_texture_features` in `src/border_texture_utils.py`*
-

Feature Extraction: Histograms

- **Pixel intensity distributions calculated for:**
 - Grayscale channel
 - Hue (H) channel
 - Saturation (S) channel
 - Value (V) channel
 - All calculated *within the masked lesion area*.
 - *Module: `calc_all_histograms` in `src/histogram_utils.py`*
-


3. Results & Visualizations

This section showcases outputs from processing sample images. (Assume images are in a local `img/` directory for this markdown)

Pipeline Stages Visualization

- Illustrates key steps in processing an image.
- Example (`ISIC_0015719.jpg`):
 - Original RGB
 - Grayscale Conversion
 - HSV (Hue Channel shown)
 - Hair Removed (Grayscale)
 - After Hair Removal (Illumination Correction typically skipped for this view)
 - Final Binary Mask (from Otsu's thresholding)

Full Pipeline Stages for: `ISIC_0015719.jpg`



`<p align="center"> Fig 1: Pipeline stages for ISIC_0015719.jpg (generated by src/full_pipeline_display.py or src/all_features_display.py) </p>`

Feature Summaries

- Textual and visual representation of extracted *intensity and color features*.
- Example (`ISIC_0015719.jpg`):

Intensity Features for ISIC_0015719.jpg
=====

Grayscale Intensity:
Mean: 186.8315
Std: 10.0121

HSV Channels:
Hue (H):
Mean: 0.0000
Std: 0.0000
Saturation (S):
Mean: 0.0000
Std: 0.0000
Value (V):
Mean: 0.0000
Std: 0.0000

Fig 2: Intensity and HSV Color Features for ISIC_0015719.jpg (generated by src/intensity_features_display.py or similar)

Channel Histograms (Masked Region)

- Visualizes pixel distribution for different channels within the segmented lesion.
- Example (ISIC_0015719 . jpg): Shows histograms for Grayscale, Hue, Saturation, and Value.

Channel Histograms (Masked Region) for: ISIC_0015719.jpg






Fig 3: Channel Histograms for ISIC_0015719.jpg (generated by src/masked_histogram_display.py or src/full_pipeline_display.py)

More Examples:

ISIC_0052212.jpg



ISIC_1132496.jpg



4. Code Structure & Key Modules

The project is organized into several Python modules within the `src/` directory:

- `color_utils.py` : Color space conversions.
- `preprocessing.py` : Hair removal, illumination correction.
- `segmentation.py` : Thresholding methods, mask application.
- `morphology.py` : Morphological operations (opening, closing).
- `histogram_utils.py` : Masked histogram calculations.
- `feature_extraction.py` : Core logic for calculating intensity, color features.
- `border_texture_utils.py` : Calculates border and texture specific features.
- `data_loader.py` : Handles image and metadata loading.
- `batch_processor.py` : Enables processing of multiple images.

Display Scripts: Numerous scripts like `full_pipeline_display.py`, `all_features_display.py`, `intensity_features_display.py`, etc., allow visualization of different pipeline stages and extracted features.

5. Challenges & Learnings

- Segmentation Accuracy:** Achieving perfect lesion segmentation is challenging due to image variability (hair, bubbles, skin lines, low contrast). Otsu's method provides a good baseline, and morphological cleanup helps significantly.
- Feature Relevance:** Understanding how different features (intensity, color, texture, border) capture various aspects of lesion appearance.
- Pipeline Integration:** Ensuring each step correctly feeds into the next and that data formats are consistent.
- Modularity:** Structuring the code into reusable functions and modules (e.g., for color conversion, preprocessing, feature calculation) was key for development and testing.

6. Conclusion & Future Work

Achievements

- Successfully implemented an end-to-end image processing pipeline for skin lesion analysis.

- Developed modules for key DIP tasks: preprocessing, segmentation, and feature extraction.
- Extracted a diverse set of features:
 - Intensity (Grayscale statistics)
 - Color (HSV channel statistics)
 - Border (Gradient-based metrics)
 - Texture (Top-hat/Bottom-hat based metrics)
- Created comprehensive visualization tools to inspect pipeline stages and results.

Future Work

- **Advanced Segmentation:** Explore more robust segmentation techniques (e.g., active contours, watershed, machine learning-based segmentation).
 - **Implement PRD Feature Set (f1-f28):**
 - Introduce the "Ida (Darkness)" channel.
 - Implement the specific Brightness, Saturation, and Darkness features (f1-f28) involving higher-order statistics (skewness, kurtosis), entropy, etc., as outlined in the project's PRD.
 - **Machine Learning Integration:** Use the extracted features to train classification models (e.g., SVM, Random Forest, Neural Networks) to predict lesion malignancy.
 - **Quantitative Evaluation:** If detailed ground truth masks become available, perform quantitative evaluation of segmentation accuracy (e.g., Dice coefficient, Jaccard index).
 - **Expanded Dataset Processing:** Utilize the `batch_processor.py` to analyze a larger subset of the dataset.
-

Thank You & Questions?
