

Skin Lesion Analysis

By Umut Celik Email: umut.celik@cix.csi.cuny.edu

CSC 741 - Digital Image Processing Midterm Project

1. Introduction

The Challenge

- Automated analysis of dermoscopic images is crucial for early skin cancer detection.
- These images present challenges: varying illumination, hair, diverse skin tones, and subtle lesion characteristics.

Project Goal

- To develop an image processing pipeline to:
 1. Segment (isolate) skin lesions from dermoscopic images.
 2. Extract a comprehensive set of visual features from the segmented lesions.
- Apply fundamental Digital Image Processing techniques to a real-world medical imaging problem.

Dataset

- **Source:** SIIM-ISIC Melanoma Classification (Kaggle).
- **Content:** Dermoscopic images (JPEG format) and associated metadata (e.g., labels).

2. Image Processing Pipeline

A multi-step approach to process images and extract meaningful data:

Image Loading & Conversion -> Preprocessing -> Color Space Transformation -> Segmentation -> Feature Extraction

Step 1: Image Loading & Initial Conversion

- **Action:** Load JPEG images.
- **Conversion:** Convert from BGR (OpenCV default) to RGB for consistency in processing and display.
 - Files involved: `src/data_loader.py`

Step 2: Preprocessing

- **Grayscale Conversion:**
 - Convert RGB images to grayscale for intensity-based operations.
 - Function: `rgb_to_grayscale` in `src/color_utils.py`
- **Hair Removal:**
 - Technique: Morphological Black-Hat filtering to identify hair structures, followed by inpainting to remove them.
 - Reduces noise and artifacts caused by hair.
 - Function: `remove_hair` in `src/preprocessing.py`
- **(Illumination Correction):**
 - Concept: Use morphological opening to estimate and correct non-uniform background illumination.
 - Status: Implemented (`correct_illumination` in `src/preprocessing.py`), but noted as potentially skipped in some display pipelines for directness to segmentation.

Step 3: Color Space Transformation

- **RGB to HSV:**
 - Convert RGB images to Hue, Saturation, Value (HSV) color space.
 - HSV is often more intuitive for color-based feature extraction:
 - **H (Hue):** Dominant color.
 - **S (Saturation):** Purity/intensity of color.
 - **V (Value):** Brightness.
 - Function: `rgb_to_hsv` in `src/color_utils.py`
- **Ida (Darkness) Channel:**
 - Calculate Ida channel defined as $\max(R, G, B) - \min(R, G, B)$ per pixel.

- Represents color spread or darkness in the image.
 - Used for darkness-based feature extraction (f19-f28).
 - *Function: `calculate_ida_channel` in `src/color_utils.py`*
-

Step 4: Lesion Segmentation

- **Goal:** Create a binary mask that accurately isolates the lesion area.
 - **Method: Otsu's Thresholding**
 - Applied to the preprocessed (hair-removed) grayscale image.
 - Automatically determines an optimal threshold value to separate lesion pixels from background.
 - *Function: `otsu_threshold` in `src/segmentation.py`*
 - **Mask Refinement:**
 - Technique: Morphological operations.
 - **Opening:** Erosion followed by dilation (removes small noise/objects).
 - **Closing:** Dilation followed by erosion (fills small holes within the lesion).
 - Improves the quality and contiguity of the mask.
 - *Functions: `opening`, `closing` in `src/morphology.py`, utilized by `apply_threshold` in `src/segmentation.py`*
-

Step 5: Feature Extraction

Extracted from the segmented lesion area (defined by the refined mask):

Feature Extraction: Intensity Features

- **From Grayscale image:**
 - Mean pixel intensity
 - Standard deviation of pixel intensities
 - Minimum pixel intensity
 - Maximum pixel intensity
 - *Module: `calculate_intensity_stats` in `src/feature_extraction.py`*
-

Feature Extraction: Color Features (HSV)

- **From HSV channels:**
 - **Hue (H):**
 - Mean Hue (using circular mean calculation for angular data)
 - Standard deviation of Hue
 - **Saturation (S):**
 - Mean Saturation
 - Standard deviation of Saturation
 - **Value (V):**
 - Mean Value (Brightness)
 - Standard deviation of Value
 - *Module: `calculate_hsv_stats` in `src/feature_extraction.py`*
-

Feature Extraction: Border Features

- **From Grayscale image & mask:**
 - Calculated using morphological gradient on the lesion border.
 - Features include:
 - Mean, Standard Deviation, and Max of border gradient values.
 - Border Irregularity (ratio of border pixels to mask area).
 - *Module: `calculate_gradient_features` in `src/border_texture_utils.py`*
-

Feature Extraction: Texture Features

- **From Grayscale image & mask:**
 - Calculated using Morphological Top-Hat (highlights bright details) and Bottom-Hat (highlights dark details) transforms.
 - Features include:
 - Mean and Standard Deviation of Top-Hat transform values.
 - Mean and Standard Deviation of Bottom-Hat transform values.
 - Texture Contrast Index (sum of Top-Hat and Bottom-Hat means).
 - *Module: `calculate_texture_features` in `src/border_texture_utils.py`*
-

Feature Extraction: Histograms

- **Pixel intensity distributions calculated for:**

- Grayscale channel
- Hue (H) channel
- Saturation (S) channel
- Value (V) channel
- Ida (Darkness) channel
- All calculated *within the masked lesion area*.
- Module: `calc_all_histograms` in `src/histogram_utils.py`

Feature Extraction: Advanced Statistical Features (f1-f28)

- **Brightness Features (f1-f9)** from Value (V) channel:
 - Basic statistics: Mean (f1), Standard Deviation (f2)
 - Higher-order moments: Skewness (f3), Kurtosis (f4)
 - Information theory: Entropy (f5)
 - Histogram analysis: Average differences (f6), Sum of 10 largest differences (f7)
 - Range ratios: High-to-mid range (f8), Mid-to-low range (f9)
 - Module: `calculate_brightness_features` in `src/feature_extraction.py`
- **Saturation Features (f10-f18)** from Saturation (S) channel:
 - Similar statistical metrics (f10-f18) applied to Saturation channel
 - Module: `calculate_saturation_features` in `src/feature_extraction.py`
- **Darkness Features (f19-f28)** from Ida channel:
 - Similar statistical metrics (f19-f27) applied to Ida (Darkness) channel
 - Additional feature: Bounding box coverage (f28)
 - Module: `calculate_darkness_features` in `src/feature_extraction.py`


3. Results & Visualizations

This section showcases outputs from processing sample images. (Assume images are in a local `img/` directory for this markdown)

Pipeline Stages Visualization

- Illustrates key steps in processing an image.
- Example (`ISIC_0015719.jpg`):
 - Original RGB
 - Grayscale Conversion
 - HSV (Hue Channel shown)
 - Hair Removed (Grayscale)
 - After Hair Removal (Illumination Correction typically skipped for this view)
 - Final Binary Mask (from Otsu's thresholding)

Full Pipeline Stages for: `ISIC_0015719.jpg`



`<p align="center"> Fig 1: Pipeline stages for ISIC_0015719.jpg (generated by src/full_pipeline_display.py or src/all_features_display.py) </p>`

Feature Summaries

- Textual and visual representation of extracted *intensity and color features*.
- Example (`ISIC_0015719.jpg`):

Intensity Features for `ISIC_0015719.jpg`

=====

Grayscale Intensity:
Mean: 186.8315
Std: 10.0121

HSV Channels:
Hue (H):
Mean: 0.0000
Std: 0.0000
Saturation (S):
Mean: 0.0000
Std: 0.0000
Value (V):
Mean: 0.0000
Std: 0.0000

`<p align="center"> Fig 2: Intensity and HSV Color Features for ISIC_0015719.jpg (generated by src/intensity_features_display.py or similar) </p>`

Channel Histograms (Masked Region)

- Visualizes pixel distribution for different channels within the segmented lesion.
- Example (`ISIC_0015719.jpg`): Shows histograms for Grayscale, Hue, Saturation, and Value.

Channel Histograms (Masked Region) for: ISIC_0015719.jpg





Fig 3: Channel Histograms for ISIC_0015719.jpg (generated by `src/masked_histogram_display.py` or `src/full_pipeline_display.py`)

Advanced Features Visualization (f1-f28)


- Visualizes the newly implemented advanced features including the Ida channel.
- Example (ISIC_0015719 . jpg): Shows Value (V), Saturation (S), and Ida channels along with segmentation and corresponding histograms.

Advanced Features Visualization (f1-f28) - ISIC_0015719




More Examples:

ISIC_0052212.jpg




Advanced Features Visualization (f1-f28) - ISIC_0052212



ISIC_1132496.jpg




Original Image




Segmentation Mask




Ida (Darkness) Channel




Value (V) Channel (f1-f9)



Saturation (S) Channel (f10-f18)



Advanced Features (f1-f28)



Advanced Features (f1-f28)

4. Code Structure & Key Modules

The project is organized into several Python modules within the `src/` directory:

- `color_utils.py`: Color space conversions, Ida channel calculation.
- `preprocessing.py`: Hair removal, illumination correction.
- `segmentation.py`: Thresholding methods, mask application.
- `morphology.py`: Morphological operations (opening, closing).
- `histogram_utils.py`: Masked histogram calculations, advanced statistical functions.
- `feature_extraction.py`: Core logic for calculating features (intensity, color, advanced f1-f28).
- `border_texture_utils.py`: Calculates border and texture specific features.
- `data_loader.py`: Handles image and metadata loading.
- `batch_processor.py`: Enables processing of multiple images.
- `advanced_features_display.py`: Visualizes advanced features (f1-f28).

Display Scripts: Numerous scripts like `full_pipeline_display.py`, `all_features_display.py`, `intensity_features_display.py`, etc., allow visualization of different pipeline stages and extracted features.

5. Challenges & Learnings

- **Segmentation Accuracy:** Achieving perfect lesion segmentation is challenging due to image variability (hair, bubbles, skin lines, low contrast). Otsu's method provides a good baseline, and morphological cleanup helps significantly.
- **Feature Relevance:** Understanding how different features (intensity, color, texture, border) capture various aspects of lesion appearance.
- **Pipeline Integration:** Ensuring each step correctly feeds into the next and that data formats are consistent.
- **Modularity:** Structuring the code into reusable functions and modules (e.g., for color conversion, preprocessing, feature calculation) was key for development and testing.

6. Conclusion & Future Work

Achievements

- Successfully implemented an end-to-end image processing pipeline for skin lesion analysis.
- Developed modules for key DIP tasks: preprocessing, segmentation, and feature extraction.
- Extracted a comprehensive set of features:

- Intensity (Grayscale statistics)
- Color (HSV channel statistics)
- Border (Gradient-based metrics)
- Texture (Top-hat/Bottom-hat based metrics)
- Advanced Statistical Features (f1-f28):
 - Brightness Features (Value channel analysis)
 - Saturation Features (Saturation channel analysis)
 - Darkness Features (Ida channel analysis)
- Created comprehensive visualization tools to inspect pipeline stages and results.
- Implemented batch processing capabilities for analyzing multiple images.

Future Work

- **Advanced Segmentation:** Explore more robust segmentation techniques (e.g., active contours, watershed, machine learning-based segmentation).
 - ✓ **Implemented PRD Feature Set (f1-f28):**
 - ✓ Introduced the "Ida (Darkness)" channel.
 - ✓ Implemented the specific Brightness (f1-f9), Saturation (f10-f18), and Darkness (f19-f28) features involving higher-order statistics (skewness, kurtosis), entropy, etc., as outlined in the project's PRD.
 - **Feature Analysis:** Evaluate the discriminative power of each feature for melanoma detection, particularly the newly implemented advanced features.
 - **Machine Learning Integration:** Use the extracted features to train classification models (e.g., SVM, Random Forest, Neural Networks) to predict lesion malignancy.
 - **Quantitative Evaluation:** If detailed ground truth masks become available, perform quantitative evaluation of segmentation accuracy (e.g., Dice coefficient, Jaccard index).
 - **Expanded Dataset Processing:** Utilize the `batch_processor.py` to analyze a larger subset of the dataset and generate comprehensive feature statistics.
-

Thank You & Questions?
