

- Skin Lesion Analysis: Segmentation & Feature Extraction
  - CSC 741 - Digital Image Processing Midterm Project
  - 1. Introduction
    - The Challenge
    - Project Goal
    - Dataset
  - 2. Image Processing Pipeline
    - Step 1: Image Loading & Initial Conversion
    - Step 2: Preprocessing
    - Step 3: Color Space Transformation
    - Step 4: Lesion Segmentation
    - Step 5: Feature Extraction
      - Feature Extraction: Intensity Features
      - Feature Extraction: Color Features (HSV)
      - Feature Extraction: Border Features
      - Feature Extraction: Texture Features
      - Feature Extraction: Histograms
  - 3. Results & Visualizations
    - Pipeline Stages Visualization
    - Feature Summaries
    - Channel Histograms (Masked Region)
    - More Examples:
  - 4. Code Structure & Key Modules
  - 5. Challenges & Learnings
  - 6. Conclusion & Future Work
    - Achievements
    - Future Work
  - Thank You & Questions?

# Skin Lesion Analysis: Segmentation & Feature Extraction

---

By **Umut Celik** Email: [umut.celik@cix.csi.cuny.edu](mailto:umut.celik@cix.csi.cuny.edu)

# CSC 741 - Digital Image Processing

## Midterm Project

---

## 1. Introduction

---

### The Challenge

- Automated analysis of dermoscopic images is crucial for early skin cancer detection.
- These images present challenges: varying illumination, hair, diverse skin tones, and subtle lesion characteristics.

### Project Goal

- To develop an image processing pipeline to:
  1. Segment (isolate) skin lesions from dermoscopic images.
  2. Extract a comprehensive set of visual features from the segmented lesions.
- Apply fundamental Digital Image Processing techniques to a real-world medical imaging problem.

### Dataset

- **Source:** SIIM-ISIC Melanoma Classification (Kaggle).
  - **Content:** Dermoscopic images (JPEG format) and associated metadata (e.g., labels).
- 

## 2. Image Processing Pipeline

---

A multi-step approach to process images and extract meaningful data:

**Image Loading & Conversion -> Preprocessing -> Color Space Transformation -> Segmentation -> Feature Extraction**

---

# Step 1: Image Loading & Initial Conversion

- **Action:** Load JPEG images.
  - **Conversion:** Convert from BGR (OpenCV default) to RGB for consistency in processing and display.
    - *Files involved:* `src/data_loader.py`
- 

## Step 2: Preprocessing

- **Grayscale Conversion:**
    - Convert RGB images to grayscale for intensity-based operations.
    - *Function:* `rgb_to_grayscale` in `src/color_utils.py`
  - **Hair Removal:**
    - Technique: Morphological Black-Hat filtering to identify hair structures, followed by inpainting to remove them.
    - Reduces noise and artifacts caused by hair.
    - *Function:* `remove_hair` in `src/preprocessing.py`
  - **(Illumination Correction):**
    - *Concept:* Use morphological opening to estimate and correct non-uniform background illumination.
    - *Status:* Implemented (`correct_illumination` in `src/preprocessing.py`), but noted as potentially skipped in some display pipelines for directness to segmentation.
- 

## Step 3: Color Space Transformation

- **RGB to HSV:**
  - Convert RGB images to Hue, Saturation, Value (HSV) color space.
  - HSV is often more intuitive for color-based feature extraction:
    - **H (Hue):** Dominant color.
    - **S (Saturation):** Purity/intensity of color.
    - **V (Value):** Brightness.
  - *Function:* `rgb_to_hsv` in `src/color_utils.py`

---

## Step 4: Lesion Segmentation

- **Goal:** Create a binary mask that accurately isolates the lesion area.
  - **Method: Otsu's Thresholding**
    - Applied to the preprocessed (hair-removed) grayscale image.
    - Automatically determines an optimal threshold value to separate lesion pixels from background.
    - *Function: `otsu_threshold` in `src/segmentation.py`*
  - **Mask Refinement:**
    - Technique: Morphological operations.
      - **Opening:** Erosion followed by dilation (removes small noise/objects).
      - **Closing:** Dilation followed by erosion (fills small holes within the lesion).
    - Improves the quality and contiguity of the mask.
    - *Functions: `opening`, `closing` in `src/morphology.py`, utilized by `apply_threshold` in `src/segmentation.py`*
- 

## Step 5: Feature Extraction

Extracted from the segmented lesion area (defined by the refined mask):

---

### Feature Extraction: Intensity Features

- **From Grayscale image:**
    - Mean pixel intensity
    - Standard deviation of pixel intensities
    - Minimum pixel intensity
    - Maximum pixel intensity
    - *Module: `calculate_intensity_stats` in `src/feature_extraction.py`*
- 

### Feature Extraction: Color Features (HSV)

- **From HSV channels:**
  - **Hue (H):**

- Mean Hue (using circular mean calculation for angular data)
  - Standard deviation of Hue
  - **Saturation (S):**
    - Mean Saturation
    - Standard deviation of Saturation
  - **Value (V):**
    - Mean Value (Brightness)
    - Standard deviation of Value
  - *Module: `calculate_hsv_stats` in `src/feature_extraction.py`*
- 

## Feature Extraction: Border Features

- **From Grayscale image & mask:**
    - Calculated using morphological gradient on the lesion border.
    - Features include:
      - Mean, Standard Deviation, and Max of border gradient values.
      - Border Irregularity (ratio of border pixels to mask area).
    - *Module: `calculate_gradient_features` in `src/border_texture_utils.py`*
- 

## Feature Extraction: Texture Features

- **From Grayscale image & mask:**
    - Calculated using Morphological Top-Hat (highlights bright details) and Bottom-Hat (highlights dark details) transforms.
    - Features include:
      - Mean and Standard Deviation of Top-Hat transform values.
      - Mean and Standard Deviation of Bottom-Hat transform values.
      - Texture Contrast Index (sum of Top-Hat and Bottom-Hat means).
    - *Module: `calculate_texture_features` in `src/border_texture_utils.py`*
- 

## Feature Extraction: Histograms

- **Pixel intensity distributions calculated for:**
  - Grayscale channel

- Hue (H) channel
  - Saturation (S) channel
  - Value (V) channel
  - All calculated *within the masked lesion area*.
  - Module: `calc_all_histograms` in `src/histogram_utils.py`
- 


## 3. Results & Visualizations

---

This section showcases outputs from processing sample images. (*Assume images are in a local `img/` directory for this markdown*)

### Pipeline Stages Visualization

- Illustrates key steps in processing an image.
- Example (`ISIC_0015719.jpg`):
  - Original RGB
  - Grayscale Conversion
  - HSV (Hue Channel shown)
  - Hair Removed (Grayscale)
  - After Hair Removal (Illumination Correction typically skipped for this view)
  - Final Binary Mask (from Otsu's thresholding)



*Fig 1: Pipeline stages for ISIC\_0015719.jpg (generated by `src/full_pipeline_display.py` or `src/all_features_display.py`)*

---

## Feature Summaries

- Textual and visual representation of extracted *intensity and color features*.
- Example ([ISIC\\_0015719.jpg](#)):

#### Intensity Features for ISIC\_0015719.jpg

##### Grayscale Intensity:

Mean: 186.8315

Std: 10.0121

##### HSV Channels:

Hue (H):

Mean: 0.0000

Std: 0.0000

Saturation (S):

Mean: 0.0000

Std: 0.0000

Value (V):

Mean: 0.0000

Std: 0.0000

*Fig 2: Intensity and HSV Color Features for ISIC\_0015719.jpg  
(generated by `src/intensity_features_display.py` or similar)*

---

## Channel Histograms (Masked Region)

- Visualizes pixel distribution for different channels within the segmented lesion.
- Example (**ISIC\_0015719.jpg**): Shows histograms for Grayscale, Hue, Saturation, and Value.



Channel Histograms (Masked Region) for: ISIC\_0015719.jpg





Fig 3: Channel Histograms for ISIC\_0015719.jpg (generated by `src/masked_histogram_display.py` or `src/full_pipeline_display.py`)

## More Examples:

### ISIC\_0052212.jpg




Pipeline Stages


Intensity Features

Channel Histograms

### ISIC\_1132496.jpg



Intensity Features for ISC\_1132496.jpg  
 Grayscale Intensity:  
 Mean: 179.5495  
 Std: 17.2239  
 HSV Channel:  
 Hue (H):  
 Mean: 0.0000  
 Std: 0.0000  
 Saturation (S):  
 Mean: 0.0000  
 Std: 0.0000  
 Value (V):  
 Mean: 0.0000  
 Std: 0.0000



Pipeline Stages

Intensity Features

Channel Histograms

## 4. Code Structure & Key Modules

The project is organized into several Python modules within the `src/` directory:

- `color_utils.py`: Color space conversions.
- `preprocessing.py`: Hair removal, illumination correction.
- `segmentation.py`: Thresholding methods, mask application.
- `morphology.py`: Morphological operations (opening, closing).
- `histogram_utils.py`: Masked histogram calculations.
- `feature_extraction.py`: Core logic for calculating intensity, color features.
- `border_texture_utils.py`: Calculates border and texture specific features.
- `data_loader.py`: Handles image and metadata loading.
- `batch_processor.py`: Enables processing of multiple images.

**Display Scripts:** Numerous scripts like `full_pipeline_display.py`, `all_features_display.py`, `intensity_features_display.py`, etc., allow visualization of different pipeline stages and extracted features.

## 5. Challenges & Learnings

- **Segmentation Accuracy:** Achieving perfect lesion segmentation is challenging due to image variability (hair, bubbles, skin lines, low contrast). Otsu's method provides a good baseline, and morphological cleanup helps significantly.
- **Feature Relevance:** Understanding how different features (intensity, color, texture, border) capture various aspects of lesion appearance.
- **Pipeline Integration:** Ensuring each step correctly feeds into the next and that data formats are consistent.

- **Modularity:** Structuring the code into reusable functions and modules (e.g., for color conversion, preprocessing, feature calculation) was key for development and testing.
- 

## 6. Conclusion & Future Work

---

### Achievements

- Successfully implemented an end-to-end image processing pipeline for skin lesion analysis.
- Developed modules for key DIP tasks: preprocessing, segmentation, and feature extraction.
- Extracted a diverse set of features:
  - Intensity (Grayscale statistics)
  - Color (HSV channel statistics)
  - Border (Gradient-based metrics)
  - Texture (Top-hat/Bottom-hat based metrics)
- Created comprehensive visualization tools to inspect pipeline stages and results.

### Future Work

- **Advanced Segmentation:** Explore more robust segmentation techniques (e.g., active contours, watershed, machine learning-based segmentation).
- **Implement PRD Feature Set (f1-f28):**
  - Introduce the "Ida (Darkness)" channel.
  - Implement the specific Brightness, Saturation, and Darkness features (f1-f28) involving higher-order statistics (skewness, kurtosis), entropy, etc., as outlined in the project's PRD.
- **Machine Learning Integration:** Use the extracted features to train classification models (e.g., SVM, Random Forest, Neural Networks) to predict lesion malignancy.
- **Quantitative Evaluation:** If detailed ground truth masks become available, perform quantitative evaluation of segmentation accuracy (e.g., Dice coefficient, Jaccard index).
- **Expanded Dataset Processing:** Utilize the `batch_processor.py` to analyze a larger subset of the dataset.

---

# Thank You & Questions?

---